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REQUIREMENTS FOR AN OPTIMAL INSTALLATION

To obtain the maximum performance and the best results, we recommend installing the balance according to the following requirements:


- a) Install it on a firm and stable surface.
- b) Place it on a stable base free of vibrations.
- c) Don't expose it to direct sunlight.
- d) Don't submit it to atmospheres with corrosive gases.
- e) Keep it in an environment free of dust and dirt.
- f) Work at temperatures between 15° and 30° C.
- g) Use the balance in environments with humidity lower than 95% HR.
- h) Keep the balance away from equipment that produces magnetic or acoustic interferences.
- i) Calibrate it in its definitive place of use and check it regularly with a calibration weight (optional), if necessary, recalibrate it.

Unpacking

Unpack the balance carefully. Inside the carton you will find the following components:

1. This operation manual (please, read it before using the balance).
2. A balance Series AH
3. A stainless steel pan
4. An AC/Adaptor of 12 V / 800 mA; don't use any other adaptor, it could damage the balance.
5. A CD-ROM with USB adaptor.
6. A USB cable
7. Five pieces of windshield.

VERY IMPORTANT: Keep the original package. It could be useful in the future in transportin the balance.

INSTALLATION

1. Remove the package from the carton.
2. Remove the pan from the cardboard.
3. Remove the AC/Adaptor (it is inside a white box).
3. Take out the balance from its protection and place it on its definitive place of use, which should be a flat, non-metallic, most stable and rigid surface.
4. Put the stainless steel pan on the pan support.
5. Level the balance by turning the four adjustable feet until the bubble level located on left frontal side remains in the centre of the circle.
6. Connect the AC/Adaptor to the balance's connector and then to an outlet of 220/240 V. 50 Hz..
7. Calibrate the balance with a proper calibration weight. This action is very important in order to obtain accurate weighing results.
8. For maximum accuracy, it is necessary to keep the balance connected for about 60 minutes before using it. However, the balance can also be used immediately with good results.

DESCRIPTION OF DISPLAY

Zero It lights when the display is at zero.

Tare It lights when a container or product have been tared.

Stable It appears when balance is stable.

Modo It lights when the "mode" key is pressed.

Setup It lights when the "Setup" key is pressed.

Cal It appears when the balance is in calibration mode.

 Battery indicator

Use  **to change** Used to navigate through the menu.

Press enter The display shows it when it is needed to confirm the function.


Press Esc The display shows it when you can return to previous menu or exit.

Weight units symbols: ct: Quilate lb: pound oz: ounce GN: grain

ozt: troy ounce

dwt: pennyweight

t1.T: tael
 pcs: pieces counting
 %: percentages
 °C: temperature

KEYBOARD DESCRIPTION


 To turn on or to turn off the scale.

MODE To select between weight, pieces-counting or percentage modes.

UNITS To select the weight unit. There are 8 weight units available.

PCS To enter the sample in pieces-counting function.

 To tare a container or product and then obtain the net weight.

 To reset the display to zero.

ENTER To confirm data or to print.

SETUP To select various menus.

ESC To return to previous menu or to exit.

 To move to the right.

 To move to the left.

 To move up.

 To move down.

WEIGHT UNIT SELECTION

Press the **UNITS** key to select the weight units, the display will view the options. There are 8 weight units available.

RESET TO ZERO

If the zero value varies during the weighing, press the $\rightarrow 0 \leftarrow$ key, the display will show zero again.

TARE

This function allows calculating the net weight of a product.

With display showing "zero", place a container on the pan, then press the  key, the balance will subtract the container's weight, put the product into the container and then read the product's weight on the display (net weight)

To return to normal weighing model, remove the container and press the  key.

OVERLOAD MESSAGE

If an object is placed on the pan and its weight exceeds the maximal capacity of the balance, the display will show the message "-----OL-----", at the same time, the balance will beep. Remove the object; otherwise it may damage the mechanism of the balance.

PIECES COUNTING FUNCTION

The pieces counting function is used to count the number of pieces deposited on the pan.

Place a container on the pan and press the  key.

Insert a certain number of pieces (10, 20, 50 or 100) into the container. This number of pieces is the initial sample, the higher number grants more accuracy in final counting.

Press the **MODE** key, select the desired sample number using the  or **UNITS** keys. To confirm the selected value, press the **PCS** key.

The balance will accept the sample and after a few minutes will indicate the pieces counting mode.

The display will show the number of pieces placed on the pan, add the rest of

pieces to obtain the total of pieces.

To count other pieces of a different weight unit, repeat the procedure.

To return to normal weighing mode, press the **MODE** key until the sign grams “g” appears on the right side.

PERCENTAGES FUNCTION (%)

Once selected the percentage mode by pressing the **MODE** key, use the keys **UNITS** or   to choose between 100.0% or 100.00%.

Place the product, considered as 100%, on the pan and press the **PCS** key. The symbol “Modo” will flash for a few seconds and then the balance will beep indicating the procedure is finished.

From this moment, the balance is ready to use the percentage function. When the sample value is less than 0.2d, the display shows “---or---%”. It means that the sample is too small, therefore, will need to add bigger samples.

Note: To exit from pieces-counting or percentages function, press the **MODE** key to return to normal weighing mode.

To return to pieces-counting or percentages function, press the **PCS** key. The balance will automatically memorize the previous data, so you can continue using both functions.


The data memory will be deleted automatically when you switch from one mode to the other.

PARAMETERS SETTING

Press the **SETUP** key to enter to parameters setting menu. Once it's pressed, the display will show the message CAL.

Use the   keys to choose the function.

The functions are:

CALIBRATION

When the display shows “CAL”, press the **ENTER** key to access to Auto-Calibration mode.

Use the keys   and   to enter the calibration weight value, then press the **ENTER** key to confirm.

Place the calibration weight on the pan, when the balance is stable, it will return to initial mode of weighing. The calibration has been finished.

FAST CALIBRATION

Keep pressed the **ENTER** key for 3 seconds. The display will show “the weight calibration value”. Put on the pan a calibration weight which corresponds to the value indicated by the balance. When the symbol “CAL” disappears, the calibration will be finished.

MAIN WEIGHT UNIT SELECTION

Press the keys **UNITS** or   to select the weight unit, then press the **ENTER** key to confirm.

To return to previous menu or to exit, press the **ESC** key.

(Default value: “g”)

SECOND UNIT WEIGHT SELECTION

(total: 8 weight units available)

Press the keys **UNITS** or   to select the weight units and press   to choose “ON” or “OFF”, then press **ENTER** to confirm.

To return to previous menu or to exit, press the **ESC** key.

AUTOPOWER OFF CONFIGURATION

Press the  key until the screen shows the message “A.OFF”, then press **ENTER**.

Use the   keys to select the auto power time, example: no, 2, 5, or 8 minutes. To confirm the desired value, press **ENTER**

To return to previous menu or to exit, press the **ESC** key.

(Default value: “no”)

DATA TRANSMISSION CONFIGURATION

Press the  key until the screen shows the message “Trn”, then press **ENTER** key.

Press   to select a “stable transmission”. Press the **ENTER** key to confirm and the **ESC** key to return to previous menu or to exit.

(Default value: “CLOSE”)
“CLOSE”:

“**stb**”: stable transmission

“**tb**”: transmission by pressing the Enter key.

“**ser**”: series transmission .

“**disp**”: transmission to an auxiliary display (optional)

BAUD RATE CONFIGURATION

Press the  key until the screen shows the message "rate", then press the **ENTER** key.

Use the keys   to select the baud rate (1200, 2400, 4800 and 9600) and press **ENTER** to confirm.

Press the **ESC** key to return to previous menu or to exit.

(Default value: 2400)

FILTER SELECTION

Press the  key until the screen shows the message "FI 105 ", then press the **ENTER** key.

Use the keys   to select the range of stability (from 1 to 15), then press **ENTER** to confirm.

Press the **ESC** key to return to previous menu or to exit.

ZERO CORRECTION CONFIGURATION

Press the  until the screen shows the message "ZEro1", then press the **ENTER** key.

Use the key   to select the zero range (from 0 to 3 divisions). To confirm, press the **ENTER** key.

Press the **ESC** key to return to previous menu or to exit.

ZEro1: correction of one division.

ZEro2: correction of two divisions.

ZEro3: correction of three divisions.

ZEro0: without correction.

AUTO-ZERO SELECTION

Press the  key until the screen shows the message "trAC1", then press the **ENTER** key

Use the   keys to select the auto-zero tracking.

Press the **ENTER** key to confirm.

Press the **ESC** key to return to previous menu or to exit.

(Default value: "2")

trAC0: it corrects deviations less than 1/3 of division.

trAC1: it corrects deviations less than 2/3 of division.

trAC2: it corrects deviations less than 1 division.

trAC3: it corrects deviations less than 4/3 of division.

trAC4: it corrects deviations less than 5/3 of division.

trAC5: it corrects deviations less than 2 divisions.

AUTO-SAMPLE SELECTION IN PIECES COUNTING

When the balance is in pieces-counting mode, this function allows making an automatic small sampling when the quantity of pieces loaded on the pan are inferior or equal to the initial quantity.

Press the  key until the screen shows the message "A.C Y", then press the **ENTER** key.

Use the   keys to select the desired mode (Yes "Y" or Not "N"). Press the **ENTER** key to confirm.

Press the **ESC** key to return to previous menu or to exit.

(Default value: "YES")

A.C Y: automatic auto-sample

A.C N: auto-sample disabled.

BACKLIGHT CONFIGURATION

Press the  key until the screen shows the message "b1.1", then press the **ENTER** key.

Use the   key to select the backlight mode, there are three options: 1, 2 and 3). Press the **ENTER** key to confirm.

Press the **ESC** key to return to previous menu or to exit.

(Default value: bl.1)

“bl.1”= Automatic backlight
 “bl.2”= Backlight always on
 “bl.3”= Backlight off

ERROR MESSAGES

The error messages are: “E1 . E2 . E5 . 0L”.
 These messages can appear if the platter is not installed properly or when the environment is not adequate. Please, contact with your distributor.

OUTPUT RS-232C

Mode: signal EIA-RS-232 C's UART, or signal USB.

Format:

Baud rate: 1200, 2400, 4800 and 9600 bauds

Data Bits 7

Parity Bit: 1


Stop Bit: 1

Code Bit: ASCII

Input Pin: 2 (RxD)

Output Pin: 3 (TxD)

Signal Pin: 5 (Masa)

CABLE TO CONNECT THE PC TO THE BALANCE

PC	BALANCE
DELTA	DELTA
9 Pins	9 Pins
Female	Male
3	2 (RxD)
2	3 (TxD)
5	5 (Grd)

DATA FORMAT

HEAD1, HEAD2, DATA UNIDAD CR

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20

HEAD1 (2 BYTES) HEAD2 (2BYTES)
OL-Overload, Mode NT-NET
ST-Stable Display GS-Gross Weight
US-Unstable Display
DATA (8 BYTES)
2D (HEX)="-" (LESS) 20 (HEX)=" "(SPACE)
2E (HEX)="." (DECIMAL POINT)
UNIT (4BYTES)
g-20 (HEX) ; 20 (HEX) ; 20 (HEX) ; 67 (HEX)
lb-20 (HEX) ; 20 (HEX) ; 6c (HEX) ; 62 (HEX)
Ti.T-74 (HEX) ; 6C (HEX) ; 2E (HEX) ; 54 (HEX)

Example of transmission

EX+0.876g, when is stable and the weight is a net weight:

HEAD, HEAD, DATA UNIDAD CR

ST , NT , + 0.876 g 0D 0A

EX-1.568lb, when is instable and the weight is a net weight.

HEAD, HEAD, DATA UNIDAD CR

US , NT , -1.568 lb 0D 0A

EX+15.24tl.T, when is stable and the weight is a net weight:

HEAD, HEAD, DATA UNIDAD CR

ST, NT , +15.24 tl.T 0D 0A

CONVERSION TABLE OF WEIGHT UNITS

1 ct (MET.CARAT) =0.1999694 g
1 lb (AVOIDUPOIS POUND) =453.59237 g
1 oz (AVOIDUPOIS OUNCE) =28.349523125g
1 GN (GRAIN) (U.K) =0.06479891 g
1 ozt (TROY OUNCE) =31.1034768 g
1 dwt (PENNYWEIGHT) =1.55517384 g
1 t1.T (TAEL) (TWN) =37.799375 g

CAPACITIES TABLE

	150g×0.005g	300g×0.01g	600g×0.01g	1500g×0.05g	3000g×0.1g
g	150.045	300.09	600.18	1500.45	3000.9

ct	750.18	1500.45	3000.9	7501.8	15004.5
lb	0.33018	0.66045	1.3009	3.3018	6.6045
oz	5.2018	10.5045	21.009	52.018	100.045
GN	2300.9	4601.8	9204.5	23009	46018
ozt	4.8018	9.6045	19.009	48.018	96.045
dwt	96.045	190.09	380.18	960.45	1900.9
tl.T	4.0018	8.0045	16.009	40.018	80.045

TECHNICAL CHARACTERISTICS

MODEL	AH-150	AH-300	AH-600	AH-1200	AH-3000
CAPACITY	150 g	300 g	600 g	1200 g	3000 g
READABILITY	0,005 g	0,01 g	0,01 g	0,02 g	0,1 g
CONSTRUCTION	Housing in ABS, pan in stainless steel.				
TARE RANGE	Full capacity				
ZERO RANGE	< + 2% maximum capacity				
WEIGHT UNITS	< + 2% maximum capacity				
FUNCTION	Weight, counting pieces and percentages				
DISPLAY	6 LCD digits (40x95 mm) with backlight				
STABILIZATION TIME	< 2 seconds				
OPERATING TEMPERATURE	0 to 40°C / 32 to 104 °F				
HUMIDITY RANGE	< 90% relative humidity, not condensation				
POWER	AC adaptor 12V DC / 500 mA and internal rechargeable battery.				
CALIBRATION	External				
BATTERY LIFE	80 hours of continuous use and 12 hours of charging time				
INTERFACE	RS-232 or USB				
PAN SIZE	116 mm (circular) 124x144 mm (square)				

GUARANTEE

This balance is guaranteed for one year from the delivery date. The guarantee covers any fabrication defect of the material.

During this period, GRAM PRECISION covers the manpower and the spare parts costs necessary for the reparation of the scale.

This guarantee does not cover the failures caused by an inappropriate use or overload.

The guarantee does not cover the freight cost (transport) necessary to repair the balance.

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